

Product datasheet for **GA118948**

Human DNLZ activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DNLZ CRISPRa kit - CRISPR gene activation of human DNL-type zinc finger
Symbol:	DNLZ
Components:	GA118948G1, DNLZ gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118948G2, DNLZ gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118948G3, DNLZ gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5SXM8</u>
Synonyms:	bA413M3.2; C9orf151; HEP; HEP1; TIMM15; ZIM17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)