

Product datasheet for **GA118953**

Human **CNTNAP3B** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CNTNAP3B CRISPRa kit - CRISPR gene activation of human contactin associated protein like 3B
Symbol:	CNTNAP3B
Components:	GA118953G1, CNTNAP3B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118953G2, CNTNAP3B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118953G3, CNTNAP3B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96NU0
Synonyms:	DKFZp686I15204
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

