

Product datasheet for **GA118964**

Human EIF3CL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3CL CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit C like
Symbol:	EIF3CL
Components:	GA118964G1, EIF3CL gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118964G2, EIF3CL gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118964G3, EIF3CL gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B5ME19</u>
Synonyms:	eukaryotic translation initiation factor 3, subunit 8, 110kDa-like; eukaryotic translation initiation factor 3, subunit C-like; OTTHUMP00000122518
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

