

Product datasheet for **GA118997**

Human TP53TG3D activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TP53TG3D CRISPRa kit – CRISPR gene activation of human TP53 target 3D
Symbol:	TP53TG3D
Components:	GA118997G1, TP53TG3D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118997G2, TP53TG3D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118997G3, TP53TG3D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9ULZ0</u>
Synonyms:	TP53TG3; TP53TG3B; TP53TG3C; TP53TG3E; TP53TG3F
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

