

Product datasheet for **GA119008**

Human **CASTOR2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CASTOR2 CRISPRa kit - CRISPR gene activation of human cytosolic arginine sensor for mTORC1 subunit 2
Symbol:	CASTOR2
Components:	GA119008G1, CASTOR2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119008G2, CASTOR2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119008G3, CASTOR2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A6NHX0</u>
Synonyms:	CASTOR2; GATSL1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

