

Product datasheet for **GA119014**

Human C6orf35 (TMEM242) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM242 CRISPRa kit - CRISPR gene activation of human transmembrane protein 242
Symbol:	TMEM242
UniProt ID:	Q9NWH2
Synonyms:	BM033; C6orf35
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)