

Product datasheet for **GA119019**

Human PRR23A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRR23A CRISPRa kit – CRISPR gene activation of human proline rich 23A
Symbol:	PRR23A
Components:	GA119019G1, PRR23A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119019G2, PRR23A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119019G3, PRR23A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A6NEV1</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

