

Product datasheet for **GA119030**

Human SHISA7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SHISA7 CRISPRa kit – CRISPR gene activation of human shisa family member 7
Symbol:	SHISA7
Components:	GA119030G1, SHISA7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119030G2, SHISA7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119030G3, SHISA7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6NL88
Synonyms:	DKFZp547H0214; FLJ34341; FLJ37640; FLJ45290; MGC120935; MGC138402
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

