

Product datasheet for **GA119037**

Human ZNF814 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF814 CRISPRa kit – CRISPR gene activation of human zinc finger protein 814
Symbol:	ZNF814
Components:	GA119037G1, ZNF814 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119037G2, ZNF814 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119037G3, ZNF814 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B7Z6K7</u>
Synonyms:	FLJ39082
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

