

Product datasheet for **GA119042**

Human TMEM229A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM229A CRISPRa kit - CRISPR gene activation of human transmembrane protein 229A
Symbol:	TMEM229A
Components:	GA119042G1, TMEM229A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119042G2, TMEM229A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119042G3, TMEM229A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B2RXF0</u>
Synonyms:	MGC189723
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)