

Product datasheet for **GA119043**

Human **TEX50** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TEX50 CRISPRa kit – CRISPR gene activation of human testis expressed 50
Symbol:	TEX50
Components:	GA119043G1, TEX50 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119043G2, TEX50 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119043G3, TEX50 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A0A1B0GTY4</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

