

Product datasheet for **GA119044**

Human LOC730183 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LOC730183 CRISPRa kit – CRISPR gene activation of human uncharacterized LOC730183
Symbol:	LOC730183
Components:	GA119044G1, LOC730183 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119044G2, LOC730183 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119044G3, LOC730183 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

