

Product datasheet for **GA119069**

Human TMEM170B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM170B CRISPRa kit – CRISPR gene activation of human transmembrane protein 170B
Symbol:	TMEM170B
Components:	GA119069G1, TMEM170B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119069G2, TMEM170B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119069G3, TMEM170B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T4T1</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

