

Product datasheet for **GA119077**

Human TRAPPC3L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRAPPC3L CRISPRa kit – CRISPR gene activation of human trafficking protein particle complex 3 like
Symbol:	TRAPPC3L
Components:	GA119077G1, TRAPPC3L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119077G2, TRAPPC3L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119077G3, TRAPPC3L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T2I5</u>
Synonyms:	bA259P20.2; BET3L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

