

## Product datasheet for **GA119083**

### Human **SIGLEC10-AS1** activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                                                 |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                                    |
| <b>Description:</b>  | LOC100129083 CRISPRa kit – CRISPR gene activation of human uncharacterized LOC100129083                                                                                                                                                                                                                                         |
| <b>Symbol:</b>       | SIGLEC10-AS1                                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>   | GA119083G1, SIGLEC10-AS1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA119083G2, SIGLEC10-AS1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA119083G3, SIGLEC10-AS1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                                         |

#### Product images:

