

Product datasheet for **GA119090**

Human XP32 (C1orf68) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C1orf68 CRISPRa kit - CRISPR gene activation of human chromosome 1 open reading frame 68
Symbol:	C1orf68
Components:	GA119090G1, XP32 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119090G2, XP32 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119090G3, XP32 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T750</u>
Synonyms:	LEP7; XP32
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

