

## Product datasheet for **GA119105**

### Human IZUMO3 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | IZUMO3 CRISPRa kit - CRISPR gene activation of human IZUMO family member 3                                                                                                                                                                                                                                    |
| <b>Symbol:</b>       | IZUMO3                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA119105G1, IZUMO3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA119105G2, IZUMO3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA119105G3, IZUMO3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>Synonyms:</b>     | bA20A20.1; C9orf134                                                                                                                                                                                                                                                                                           |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

