

Product datasheet for **GA119115**

Human **SUPT20HL1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SUPT20HL1 CRISPRa kit - CRISPR gene activation of human SUPT20H like 1
Symbol:	SUPT20HL1
Components:	GA119115G1, SUPT20HL1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119115G2, SUPT20HL1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119115G3, SUPT20HL1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	FAM48B1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

