

Product datasheet for **GA119116**

Human C17orf107 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C17orf107 CRISPRa kit - CRISPR gene activation of human chromosome 17 open reading frame 107
Symbol:	C17orf107
Components:	GA119116G1, C17orf107 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119116G2, C17orf107 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119116G3, C17orf107 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZR85</u>
Synonyms:	FLJ99427
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

