

Product datasheet for **GA119120**

Human LOC100130370 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LOC100130370 CRISPRa kit – CRISPR gene activation of human uncharacterized LOC100130370
Symbol:	LOC100130370
Components:	GA119120G1, LOC100130370 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119120G2, LOC100130370 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119120G3, LOC100130370 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	DKFZp761P0212
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

