

Product datasheet for **GA119122**

Human CD300H activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CD300H CRISPRa kit – CRISPR gene activation of human CD300H molecule (gene/pseudogene)
Symbol:	CD300H
Components:	GA119122G1, CD300H gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119122G2, CD300H gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119122G3, CD300H gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	CD300H; CD300Hs
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

