

Product datasheet for **GA119126**

Human **LRRC70** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC70 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 70
Symbol:	LRRC70
Components:	GA119126G1, LRRC70 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119126G2, LRRC70 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119126G3, LRRC70 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z2Q7</u>
Synonyms:	SLRN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)