

## Product datasheet for **GA119130**

### Human SBK3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SBK3 CRISPRa kit - CRISPR gene activation of human SH3 domain binding kinase family member 3
Symbol:	SBK3
Components:	GA119130G1, SBK3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119130G2, SBK3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119130G3, SBK3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0C264</u>
Synonyms:	SGK110
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)