

Product datasheet for **GA119131**

Human TSTD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSTD3 CRISPRa kit - CRISPR gene activation of human thiosulfate sulfurtransferase like domain containing 3
Symbol:	TSTD3
Components:	GA119131G1, TSTD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119131G2, TSTD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119131G3, TSTD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	H0UI37
Synonyms:	RP1-199J3.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)