

## Product datasheet for **GA119145**

### Human CD300LD activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | CD300LD CRISPRa kit - CRISPR gene activation of human CD300 molecule like family member d                                                                                                                                                                                                                        |
| Symbol:       | CD300LD                                                                                                                                                                                                                                                                                                          |
| Components:   | GA119145G1, CD300LD gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA119145G2, CD300LD gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA119145G3, CD300LD gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q6UXZ3</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | CD300D; CLM-4; CMRF35-A4; CMRF35A4                                                                                                                                                                                                                                                                               |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)