

Product datasheet for **GA119168**

Human DPH3P1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DPH3P1 CRISPRa kit - CRISPR gene activation of human diphthamide biosynthesis 3 pseudogene 1
Symbol:	DPH3P1
Components:	GA119168G1, DPH3P1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119168G2, DPH3P1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119168G3, DPH3P1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H4G8
Synonyms:	C20orf143; DPH3B; ZCSL1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)