

Product datasheet for **GA119169**

Human SHISAL2B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SHISAL2B CRISPRa kit - CRISPR gene activation of human shisa like 2B
Symbol:	SHISAL2B
Components:	GA119169G1, SHISAL2B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119169G2, SHISAL2B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119169G3, SHISAL2B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6NKG6
Synonyms:	family with sequence similarity 159, member B; RGD1564672
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)