

Product datasheet for **GA119181**

Human TEN1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEN1 CRISPRa kit - CRISPR gene activation of human TEN1 subunit of CST complex
Symbol:	TEN1
Components:	GA119181G1, TEN1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119181G2, TEN1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119181G3, TEN1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q86WV5
Synonyms:	C17orf106
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)