

Product datasheet for **GA119238**

Human TMEM239 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM239 CRISPRa kit - CRISPR gene activation of human transmembrane protein 239
Symbol:	TMEM239
Components:	GA119238G1, TMEM239 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119238G2, TMEM239 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119238G3, TMEM239 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WW34
Synonyms:	FLJ26142; FLJ44045
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)