

Product datasheet for **GA119253**

Human PET117 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PET117 CRISPRa kit - CRISPR gene activation of human PET117 cytochrome c oxidase chaperone
Symbol:	PET117
Components:	GA119253G1, PET117 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119253G2, PET117 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119253G3, PET117 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UWS5</u>
Synonyms:	CSRP2BP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)