

## Product datasheet for **GA119254**

### Human SPDYE2B activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | SPDYE2B CRISPRa kit - CRISPR gene activation of human speedy/RINGO cell cycle regulator family member E2B                                                                                                                                                                                                        |
| Symbol:       | SPDYE2B                                                                                                                                                                                                                                                                                                          |
| Components:   | GA119254G1, SPDYE2B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA119254G2, SPDYE2B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA119254G3, SPDYE2B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>A6NHP3</u>                                                                                                                                                                                                                                                                                                    |
| Synonyms:     | SPDYE2L                                                                                                                                                                                                                                                                                                          |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)