

Product datasheet for **GA119278**

Human **LRRC3C** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC3C CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 3C
Symbol:	LRRC3C
Components:	GA119278G1, LRRC3C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119278G2, LRRC3C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119278G3, LRRC3C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6NIW4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)