

Product datasheet for **GA119285**

Human LOC100505841 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LOC100505841 CRISPRa kit - CRISPR gene activation of human zinc finger protein 474-like
Symbol:	LOC100505841
Components:	GA119285G1, LOC100505841 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119285G2, LOC100505841 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119285G3, LOC100505841 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)