

Product datasheet for **GA119300**

Human LOC100506422 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LOC100506422 CRISPRa kit - CRISPR gene activation of human putative deoxyuridine 5'-triphosphate nucleotidohydrolase-like protein FLJ16323
Symbol:	LOC100506422
Components:	GA119300G1, LOC100506422 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119300G2, LOC100506422 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119300G3, LOC100506422 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZN92</u>
Synonyms:	FLJ16323
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)