

Product datasheet for **GA119304**

Human C16orf95 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C16orf95 CRISPRa kit - CRISPR gene activation of human chromosome 16 open reading frame 95
Symbol:	C16orf95
Components:	GA119304G1, C16orf95 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119304G2, C16orf95 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119304G3, C16orf95 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H693
Synonyms:	DKFZp434J1815; FLJ22477
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)