


Product datasheet for **GA119321**

Human TMEM178B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM178B CRISPRa kit - CRISPR gene activation of human transmembrane protein 178B
Symbol:	TMEM178B
Components:	GA119321G1, TMEM178B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119321G2, TMEM178B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119321G3, TMEM178B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	H3BS89
Synonyms:	DKFZp547G036; FLJ39162; FLJ42883
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)