


Product datasheet for **GA119399**

Human CEP295NL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CEP295NL CRISPRa kit - CRISPR gene activation of human CEP295 N-terminal like
Symbol:	CEP295NL
Components:	GA119399G1, CEP295NL gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119399G2, CEP295NL gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119399G3, CEP295NL gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96MC4
Synonyms:	DDC8; KIAA1731NL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)