

## Product datasheet for **GA119435**

### Human PYURF activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PYURF CRISPRa kit - CRISPR gene activation of human PIGY upstream reading frame
Symbol:	PYURF
UniProt ID:	<a href="#">Q96I23</a>
Synonyms:	PREY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)