

Product datasheet for **GA119448**

Human OPN1MW3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	OPN1MW3 CRISPRa kit - CRISPR gene activation of human opsin 1, medium wave sensitive 3
Symbol:	OPN1MW3
Components:	GA119448G1, OPN1MW3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA119448G2, OPN1MW3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA119448G3, OPN1MW3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0DN77</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)