

Product datasheet for **GA200228**

Mouse Birc3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Birc3 CRISPRa kit - CRISPR gene activation of mouse baculoviral IAP repeat-containing 3 |
| Symbol: | Birc3 |
| UniProt ID: | <u>O08863</u> |
| Synonyms: | Api1; Api2; AW107670; Birc2; C-IAP2; cIAP-1; cIAP-2; cIAP1; cIAP2; HIAP2; IAP1; IAP2; MIAP1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)