

Product datasheet for **GA200270**

Mouse Arg1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Arg1 CRISPRa kit - CRISPR gene activation of mouse arginase, liver
Symbol:	Arg1
UniProt ID:	Q61176
Synonyms:	AI; AI256583; Arg-1; PGIF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)