

Product datasheet for **GA200392**

Mouse **Opn1sw** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Opn1sw CRISPRa kit - CRISPR gene activation of mouse opsin 1 (cone pigments), short-wave-sensitive (color blindness, tritan)
Symbol:	Opn1sw
Components:	GA200392G1, Opn1sw gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA200392G2, Opn1sw gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA200392G3, Opn1sw gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P51491</u>
Synonyms:	AW551857; Bcp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)