

Product datasheet for **GA200475**

Mouse *Cacna1a* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Cacna1a</i> CRISPRa kit - CRISPR gene activation of mouse calcium channel, voltage-dependent, P/Q type, alpha 1A subunit
Symbol:	<i>Cacna1a</i>
UniProt ID:	P97445
Synonyms:	alpha1A; APCA; BI; <i>Caca1a</i> ; <i>Cacn1a4</i> ; Cav2.1; <i>Ccha1a</i> ; EA2; FHM; HPCA; Ia; MHP; MHP1; nmf352
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)