

Product datasheet for **GA200477**

Mouse *Cacna1c* activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | <i>Cacna1c</i> CRISPRa kit - CRISPR gene activation of mouse calcium channel, voltage-dependent, L type, alpha 1C subunit |
| Symbol: | <i>Cacna1c</i> |
| UniProt ID: | Q01815 |
| Synonyms: | Cav1.2; Cchl1a1; D930026N18Rik; MBC; MELC-CC |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)