

Product datasheet for **GA200646**

Mouse Cd9 activation kit by CRISPRa

Product data:

| | |
|---------------|-------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | Cd9 CRISPRa kit - CRISPR gene activation of mouse CD9 antigen |
| Symbol: | Cd9 |
| UniProt ID: | P40240 |
| Synonyms: | Tspan29 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)