


Product datasheet for **GA200673**

Mouse Cdk5r2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cdk5r2 CRISPRa kit - CRISPR gene activation of mouse cyclin-dependent kinase 5, regulatory subunit 2 (p39)
Symbol:	Cdk5r2
Components:	GA200673G1, Cdk5r2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA200673G2, Cdk5r2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA200673G3, Cdk5r2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q35926</u>
Synonyms:	B230310J22Rik; p39
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)