

Product datasheet for **GA200726**

Mouse Chrm4 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Chrm4 CRISPRa kit - CRISPR gene activation of mouse cholinergic receptor, muscarinic 4 |
| Symbol: | Chrm4 |
| UniProt ID: | P32211 |
| Synonyms: | Chrm-4; M4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)