

Product datasheet for **GA200768**

Mouse *Cxcr3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cxcr3 CRISPRa kit - CRISPR gene activation of mouse chemokine (C-X-C motif) receptor 3
Symbol:	Cxcr3
UniProt ID:	<u>O88410</u>
Synonyms:	Cd183; Cmkar3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)