

Product datasheet for **GA200769**

Mouse *Cxcr4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cxcr4 CRISPRa kit - CRISPR gene activation of mouse chemokine (C-X-C motif) receptor 4
Symbol:	Cxcr4
UniProt ID:	P70658
Synonyms:	b2b220Clo; CD184; Cmkar4; LESTR; PB-CKR; PBSF/SDF-1; Sdf1r
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)