

Product datasheet for **GA200773**

Mouse *Ccr3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Ccr3</i> CRISPRa kit - CRISPR gene activation of mouse chemokine (C-C motif) receptor 3
Symbol:	<i>Ccr3</i>
UniProt ID:	P51678
Synonyms:	CC-CKR3; CKR3; Cmkbr1l2; Cmkbr3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)